REF 07840

ARTHROFIX® AIR

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SIDE EFFECTS

If used appropriately and correctly, there are no general side effects, intolerances with long-term effects or allergic reactions known to date.

INTENDED PATIENT TARGET GROUP

The target group is all patients, taking into account indications and contraindications. Instruction and allocation of the correct size is carried out by medical professionals.

PUTTING IT ON AND USE

- We recommend wearing a stocking under the brace.
- Place your heel on the foot pad of the opened brace. The padded sides point upwards **1**.
- Pass the lower and upper Velcro straps through the deflection loops. Adjust the foot pad to your ankle width using the two Velcro points. Velcro the foot pad to the side shells. Now align the brace on the lower leg so that the curves of the anatomically shaped side shells with the air cushions lie over the inner or outer ankle. Hold this position **2**.
- Fix the lower Velcro fastener first. When closing the Velcro, make sure that the pre-fixed position of the side shells from step **2** is maintained **3**.
- Finally, fix the brace with the upper Velcro strap **4**.
- The stability effect of the brace is optimised by wearing it in a closed shoe.
- Die SPORLASTIC GmbH ist einem Rücknahmesystem für Verpackungen angeschlossen und kommt somit den in der Bundesrepublik Deutschland geltenden verpackungsrechtlichen Bestimmungen vollständig nach.

IMPORTANT INFORMATION

- Essentially, the indication for and duration of wearing of any orthopaedic device in general, as well as the presence of any of the conditions below in particular, should be discussed with the doctor who is in charge of your treatment.
- Check the functioning of the product with your doctor or specialist retailer.

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INSTRUCTIONS FOR USE EN

Dear patient,

We ask you to carefully follow the enclosed instructions for use. If you have any questions, please contact the doctor looking after you, your nearest specialist store or us directly.

INTENDED USE

ARTHROFIX® AIR is an ankle brace for stabilisation and / or early functional mobilisation of the ankle joint, e.g. after ankle sprains, acute capsule ruptures or chronic instabilities of the ankle joint.

INDICATIONS

- Conservative / postoperative
- All conditions requiring stabilisation of the ankle with restriction of pronation and / or supination, such as:
 - Protection following ligament suturing / reconstruction
 - Acute capsular ligament rupture
 - Functional treatment of ligament ruptures
 - Ankle fractures (Weber A)
 - Ankle sprains
 - Chronic capsular ligament instability
 - Primary and recurrence prophylaxis

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CONTRA-INDICATIONS

- Allergic, inflammatory or injury-related skin changes (e.g. swelling, redness) of the body areas to be supplied
- Impaired circulation or lymphatic soft tissue swelling
- Neurogenic disorders of the sensory system and skin trophism in the body area to be treated (sensory disorders with or without skin damage)

QUALITY MANAGEMENT SYSTEM

All products from SPORLASTIC GmbH are subject to product inspection as part of our quality management system. If you still have any complaints about our product, we kindly ask that you contact your specialist retailer.

ONLINE VERSION OF THE INSTRUCTIONS FOR USE

An online version of the instructions for use can be found on our website at www.sporlastic.de.

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