<b>REF</b>	DE Artikelnummer EN Article number FR Référence IT Numero articolo PL Numer artykułu DA Varenummer
<b>LOT</b>	DE Fertigungslosnummer, Charge EN Production batch number, batch FR Numéro de lot de fabrication IT Numero di lotto di produzione, partita PL Numer partii produkcyjnej, szarża DA Batchnummer, batch
<b>PZN</b>	DE Pharmazentralnummer EN Pharmaceutical registration number FR Numéro central de la pharmacie IT PZN (Germania) PL Numer centralny farmacji DA Farmacentralnummer

SPORLASTIC ORTHOPAEDICS

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# INSTRUCTIONS FOR USE EN

## INTENDED USE

GENUDYN® CI STEPTHRU is a frame brace for guiding and stabilising the knee joint according to the 4-point stabilisation principle with extension/flexion limitation.

## INDICATIONS

- Injuries of the knee with involvement of the ligamentous apparatus
- After ligament injuries/reconstructions ACL (anterior cruciate ligament), PCL (posterior cruciate ligament) or combined (CI)
- Instabilities with corresponding drawer signs

## CONTRA-INDICATIONS

- Allergic, inflammatory or injury-related skin changes (e.g. swelling, redness) of the body areas to be cared for.
- Impaired circulation or lymphatic soft tissue swelling.
- Neurogenic disorders of the sensory system and skin trophism in the body area to be cared for (sensory disorders with and without skin damage).

## SIDE-EFFECTS

If used appropriately and correctly, there are no general side effects, intolerances with long-term effects or allergic reactions known to date.

## INTENDED PATIENT TARGET GROUP

The target group is all patients, taking into account indications and contraindications. Instruction and allocation of the correct size is carried out by medical professionals.

## PUTTING IT ON AND USE

- **Open the click locks:** Release the four straps by pressing the click lock button on the inside of the brace. Repeat this for all four straps.
- **Putting on the brace:** Sit on the edge of a chair and place the brace on the leg with the knee bent at about 30°. The centre of the joint should be at the level of the upper edge of the kneecap (patella) 1.
- **Putting on the straps:** Pull the strap 1 tight and make sure that the strap is not twisted. Guide the strap around your lower leg. Close it by inserting it into the click closure on the inside of the brace. A short clicking sound confirms correct closure. Make sure that the strap 1 runs above the calf belly 2.
- **Tensioning the straps:** You get the necessary strap tension by loosening the strap end at the Velcro fastener and then tightening the strap so that the brace fits tightly. Then repeat these steps for straps 2 - 4. Note that when tightening strap 2, the calf cuff rises about 1 cm above the shin. This is normal and necessary for a good fit of the brace 3.
- **Closing the front Velcro fasteners:** The front calf and thigh straps only need to be adjusted by the orthopaedic technician during the first fitting 4.
- **Check:** Check the position of the brace while standing with the leg extended. The joint centre of the brace should be on the upper edge of the kneecap (patella). You can readjust the straps as needed to ensure a tight fit 5.
- **Removing the brace:** Press the click lock buttons on the inside of the brace one after the other. You can then pull the brace forwards over your foot 6.

## NOTES FOR THE ORTHOPAEDIC TECHNICIAN

- GENUDYN® CI STEPTHRU is delivered with a correct belt sequence for generating dynamic pre-tensioning at front drawer (ACL). The belt sequence must not be interchanged so that the correct closing sequence corresponds to the numbering.
- The fit, joint position and strapping must always be checked for correctness even when the patient is standing. With a correctly fitted and well adjusted brace, it is not necessary to open the Velcro fasteners of the straps. The front straps remain closed when putting on and taking off and the rear straps are opened and closed via the quick-release fasteners. The front thigh strap is not necessary for ACL-only fitting and could be removed. All straps should be closed tightly but should not cause strangulation.
- The brace frame is made of aluminium and can be individually bent by the orthopaedic technician. The orthotic joints should be as close as possible to the condyles. Due to the Living Hinge, parallel alignment of the orthotic joints is not necessary. The support on the condyles can also be individually adapted by exchanging condylar pads of different thicknesses.
- Seat belt pads are removable, should be placed centrally and can be shortened if necessary.
- All strap lengths are set during fitting and should only be changed later by the patient in case of problems, for example due to circumference changes.

### A) ADJUSTMENT IN CASE OF INVOLVEMENT OF THE ANTERIOR CRUCIATE LIGAMENT

- The numbering of the straps corresponds to the order in which they are closed. Open all straps and place the brace on the leg at about 30° flexion. The centre of the joint should be at the level of the upper edge of the patella 1A.
- Hook strap 1 into the snap lock. Adjust the strap length using the Velcro fastener on the other side of the strap. The pads can be removed and shortened if necessary. Make sure the pads are positioned centrally. Repeat this step with belt number 2; 2A.
- Now hook strap 3 into the snap lock and pull it tight so that the brace joint is pulled back about 1.5 cm behind the centre of the condyles 3A.
- Then hook strap 4 and pull it tight enough to pull the joints back to the centre of the knee 4A.
- Then close the strap at the calf. Please note that this has no digit. Finally, close the strap at the thigh (also without a digit). The straps should be tight but not constricting 5A.

### B) ADJUSTMENT IN CASE OF INVOLVEMENT OF THE POSTERIOR CRUCIATE LIGAMENT

- NOTE: The numbering of the straps corresponds to the order in which they are closed. Open all straps and place the brace on the leg at approximately 30° flexion. The centre of the joint should be at the level of the upper edge of the patella 1B.
- Hook strap 1 into the snap lock. Adjust the strap length using the Velcro fastener on the other side of the strap. The pads can be removed and shortened if necessary. Make sure the pads are positioned centrally. Repeat this step with belt number 2; 2B.
- Now the brace is adjusted once by the two unnumbered straps in such a way that a dynamic pre-tension optimally stabilises a rear drawer. Close and tighten the rear calf strap in such a way that the orthotic joints come to rest about 1.5 cm behind the centre of the condyles. With the front strap on the thigh, the orthotic joints are pulled forward again so that the joint centres are centred on the condyles 3B.
- Now hook strap 3 into the snap lock and pull it tight 4B.
- Then hook in and tighten strap 4. With the patient standing, check that the fit, joint position and belting are correct 5B.

### C) FURTHER ADJUSTMENT AND CUSTOMISATION OPTIONS:

- **Condylar pads:** The medial and lateral condylar pads should fit the knee well. In case of gaps or different fit of the medial or lateral side, use the additional condylar pads in 5 or 15 mm thickness for a perfect fit on both sides (10 mm pads are in the brace).
- The unnumbered straps are only adjusted once by the technician during fitting and should only be changed later if the circumference of the leg changes.
- **Adjusting the pad length:** For optimal fit and wearing comfort, the belt pads must be centred on the leg when the straps are tightened. If a strap pad is too long, remove it from the strap and cut it to fit 1C with scissors.
- **Adjusting the strap length:** To adjust the strap length, open the quick release. Loosen the Velcro on the other end of the strap and adjust the length accordingly 2C.
- **Adjusting the joints:** The extension wedge set to 0° is already pre-assembled when delivered. Additional wedges for extension (EXT) and flexion (FLEX) limitation are included in the delivery. In addition, the brace comes with a pin that simplifies changing the wedges. When changing the extension or flexion limitation, always check that both the medial and

lateral joint stops are set identically.

- **The following joint stops are possible:**

EXT: 0°, 10°, 20°, 30°

FLEX: 0°, 10°, 20°, 30°, 45°, 60°, 90°

Immobilisation at 0°, 10°, 20°, 30° possible

- **Adjusting the FLEX and EXT stops:** When delivered, the 0° extension wedges are already pre-assembled.
- **The wedges for limiting movement can be changed quickly and easily without tools - even without taking off the brace.**
  1. To do this, simply pull the corresponding wedges, each at both joints, out of joint **3C**.
  2. Then insert the desired wedges in the same way. The correct engagement of the wedges can be clearly felt and heard.
- **Immobilisation:** If immobilisation is desired, the wedges must be inserted in the following order (synchronously at both joints):
  1. insert FLEX wedges
  2. insert EXT wedges

Disassembly is done in reverse order.

- A small lever tool is supplied to make it easier to remove the wedges. Insert the tip into the corresponding holes in the joints and use it to lever the wedges out **4C**.

### **IMPORTANT INFORMATION**

- Essentially, the indication for and duration of wearing of any orthopaedic device in general, as well as the presence of any of the conditions below in particular, should be discussed with the doctor who is in charge of your treatment.
- Check the functioning of the product with your doctor or specialist retailer.
- The product is designed to treat one patient.
- Do not wear the product directly against open wounds.
- In order for the product to have a long service life and function, it may not be worn in conjunction with oily/greasy or acidic products, ointments or lotions.
- The product can come into contact with splash water.
- Please always close any Velcro fasteners present properly. Damage caused by improper closing of the Velcro fasteners does not constitute a reason for complaint.
- In the event of serious incidents when using this aid, immediately contact your doctor, your health care supply retailer or contact SPORLASTIC directly at [info@sporlastic.de](mailto:info@sporlastic.de) to report the incident.
- Any improper modifications to the product and/or improper use of the above-mentioned product exclude any product liability on the part of the manufacturer.
- Product may contain metallic objects. We recommend not putting them in contact with electromagnetic interference.
- Possible health-related mutual risks or other disadvantages with certain treatments that may arise in association with the use of the product must be discussed with the doctor in charge of treatment.

Please return the packaging and the product to the local recycling collection system. Please comply with local regulations.

SPORLASTIC GmbH is affiliated with a returns system for packaging and therefore fully complies with the packaging legislation in the Federal Republic of Germany.

### **CARE**

#### **Cleaning the brace**

The pads need to be removed and cleaned regularly depending on the level of activity and perspiration. All pads and cushions are held in place by Velcro fasteners and are therefore easy to remove from the frame, hinges and straps.

Wipe the pads or wash the pads by hand in a mild, antibacterial soapy water (30°C), then rinse well and allow to air dry. Do not dry the product in the tumble dryer, on the heater or in the microwave. Reattach the pads and cushions by pressing them firmly against the corresponding points.

#### **Cleaning the joints**

The joints of the brace are pre-lubricated. However, if sand, dirt or water gets into the hinges, they may need to be cleaned and lubricated.

Remove the condylar pad, loosen the hinge screws with the Allen key, clean the hinge thoroughly and then lubricate it with a commercially available synthetic lubricant. Now check that the hinge works properly. Then refit the condylar pad and the hinge screws.

### **TRANSPORT AND STORAGE CONDITIONS**

Please ensure that the product is stored in a dry place and is protected against moisture and sunlight. Store the product at a normal temperature and humidity.

### **QUALITY MANAGEMENT SYSTEM**

All products from SPORLASTIC GmbH are subject to product inspection as part of our quality management system. If you still have any complaints about our product, we kindly ask that you contact your specialist retailer.

An online version of the instructions for use can be found on our website at [www.sporlastic.de](http://www.sporlastic.de).

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