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
INSTRUCTIONS FOR USE EN

INTENDED USE

TRIAC II is a dynamic scoliosis brace with an open frame construction.

INDICATIONS

- Conservative
- Idiopathic scoliosis up to 35° Cobb

CONTRA-INDICATIONS

- Scoliosis vertex higher than Th 7
- Scoliosis vertex Th 11/12, L 4/5
- Paralysis scoliosis
- Goal: Prevention of progression of the disease through early use of braces

SIDE-EFFECTS

If used appropriately and correctly, there are no general side effects, intolerances with long-term effects or allergic reactions known to date.

INTENDED PATIENT TARGET GROUP

The target group is all patients, taking into account indications and contraindications. Instruction and allocation of the correct size is carried out by medical professionals.

PUTTING IT ON AND USE


- We recommend wearing a thin, tight-fitting and preferably seamless T-shirt under the product.
- The TRIAC II brace should be worn 23.5 hours a day in consultation with the doctor in charge and should only be taken off for personal hygiene.
- The brace may be taken off during sports activities. Sport is considered as wearing time.

⚠️ The TRIAC II brace requires an individual and qualified initial fitting to the patient by appropriately trained specialist personnel. This is the only way to guarantee performance and wearing comfort.


⚠️ We recommend practising putting on and taking off together during the fitting by the orthopaedic technician.

- To put on the brace, open the pelvic belt and the chest closure of the TRIAC II brace.
- The two elastic shoulder straps determine the position of the TRIAC II brace.
- The arm opposite the open chest closure is the first to dip through the loop of the shoulder strap above the joint pocket.
- Then the head dives through the opening between the two shoulder straps.
- The TRIAC II brace is now open at the side of the body. The two shoulder straps are in position and the crescent-shaped chest pads are at the front without tension.
- Next, the chest closure is closed. To do this, hook the fastener with the hole onto the crescent-shaped chest pad.
- Finally, the pelvic belt is closed under tension. The correct tension as determined by the orthopaedic technician is created when the metal loop of the pelvic belt comes to lie exactly between the fleece and Velcro straps of the belt.
- After the TRIAC II brace has been put on correctly, the T-shirt underneath should be stretched without wrinkles to prevent possible pressure points.
- The TRIAC II brace is removed in reverse order.
- Especially during the first weeks of therapy with the TRIAC II brace, the posture changes, which is why a fit check is necessary at shorter intervals. If necessary, strap tensions have to be corrected and thus the corrective effect of the TRIAC II brace has to be adapted to the respective stage of development.
- Please arrange follow-up appointments with your orthopaedic technician during the initial fitting. A service pass is included in the scope of delivery.
- It is advisable to arrange a check-up appointment with the prescribing doctor after about two weeks of getting used to the TRIAC II brace, who should check the fit and effect of the TRIAC II brace.

INFORMATION FOR ORTHOPAEDIC TECHNICIANS

- | | |
|--|-------------------------------------|
| A Plastic shell with pad button | E Elastic strap |
| B Plastic shell | E1 Elastic strap |
| C TRIAC II - joint bag | F Thoracic pad |
| C1 Belt | F1 Connecting rod lumbar pad |
| C2 Belt | G Thoracic pad |
| D Upper belts | G1 Velcro tape |

- | | |
|----------------------------------|----------------------------|
| H Lumbar pad | I1 Velcro fasteners |
| H1 Connecting rod | J Lumbar cushion |
| H2 Lumbar lever | J1 Velcro tape |
| H3 Cushion for J3 | J2 Cushion for H2 |
| H4 Plastic flap | J3 Lumbar lever |
| I Velcro fasteners lumbar | K Bracket |

 You can request a dimension sheet for the product from us or download it directly from our homepage at www.sporlastic.de.

• FIGURE 1-3

Heat plastic shells **A** and **B** in a hot water bath to approx. 100°C until they become flexible, then mould them laterally to the ribs as shown in Fig. 1 + 2. Please make sure that there is enough space under the armpit for the joint pocket **C** in the plastic shell **B**.

• FIGURE 4-5

Measurements are taken from hole to hole, the upper and lower holes respectively, as shown in Fig. 6. Make sure that the plastic shells **A** and **B** do not slip when measuring. The measured dimensions, as shown in picture 4, are then carefully transferred to the straps **E** + **E1**. Thread the elastic straps **E** + **E1** into the upper slots of the connecting rod **F1**.

• FIGURE 6

Then, as shown in Fig. 6, rivet the straps **C1** + **C2** to the plastic shells **A** and **B**. Make sure that the elastic support straps **E** + **E1** are riveted under the upper strap **D**.

• FIGURE 7

Adjust the straps on the patient in such a way that the connecting rod **F1** does not run centrally but to the side of the spine.

• FIGURE 8

When all straps are correctly positioned, the force on the thoracic pad **F** is increased. Since the tension of the thoracic pelotte pad **F** is adjusted in the open state, it is advisable to unhook the thoracic pelotte pad **F** from the pelotte head of the plastic shell **A** before adjusting it and only re-hook it after adjustment. After this procedure, the scale on the thoracic pelotte pad should be at approx. 2 to 2.5 in the exhaled state.

• FIGURE 9

Relief cuts are made on the lower edge of the lumbar pad **H** with the hole punch and the scissors. This reduces the distal edge pressure of the lumbar pad **H**. By setting the connecting rod **H1** outwards (approx. 30° to 45°), forces are built up on the lumbar and thoracic pelotte. The more the connecting rod **H1** is bent outwards, the more corrective force the TRIAC II scoliosis brace exerts on the spine.

• FIGURE 10

The orthopaedic technician must attach a custom-made lumbar pad **J**, e.g. by means of Velcro tape **J1**, to the lumbar support. Pad and Velcro are not included with the product.

• FIGURE 11

Cut the connecting rod **H1** to the desired length. Make sure that the rod is not too short. It should reach at least to the upper edge of the joint pocket **C**. It is important that the patient looks straight ahead and stands upright to avoid cutting the connecting rod **H1** too short.

• FIGURE 12-15

The lumbar levers **H2** and **J3** are executed dorsally into the middle of the gluteus maximus muscle and pulled from ventrally to the middle of the thigh. Velcro straps **I** and **I1** are then attached with clamps so that they run horizontally above the trochanter. During the sitting test, the ventral lumbar lever is cut to length so that the patient can sit well and bend forward without obstruction.

IMPORTANT: The lumbar levers **H2** and **J3** should not be cut unnecessarily short to ensure maximum leverage.

• FIGURE 16

The metal connections of the lumbar pelotte pad **H** are restricted in such a way that pressure points are avoided.

• FIGURE 17

The more the connecting rod **H1** is bent outwards, the more correction the TRIAC II scoliosis brace exerts on the spine.

• FIGURE 18

There should be at least 1 finger width of space between the joint pocket **C** and the patient to ensure sufficient correction clearance. When sufficient clearance can be felt at the joint pocket and the scale of the thoracic support **F** is at 2.5 - 3, it can be riveted to the connecting rod of the lumbar support **F1**.

• FIGURE 19

Use the soldering iron to melt holes into the measured straps and secure them with hollow rivets. The elastic beams **E** + **E1** should not be secured with hollow rivets to allow for later readjustment.

• FIGURE 20-21

Finally, cut the plastic tab **H4** to length and drill the connecting rod **H1** at the upper end. Then rivet them both together and thread them into the tab of the hinge pocket **C**. This is then secured with a hollow rivet so that the connecting rod **H1** can move freely in the joint pocket **C**.

• FIGURE 22

If required, a custom-made pad **G** can be attached to the thoracic pad **F**, e.g. by means of Velcro **G1**. Pad and Velcro are not included with the product.

• FIGURE 23

To put on the TRIAC II scoliosis brace, first close the thoracic closure, then the lumbar closure. To remove the TRIAC II scoliosis brace, proceed in reverse order.

⚠ IMPORTANT INFORMATION

- Essentially, the indication for and duration of wearing of any orthopaedic device in general, as well as the presence of any of the conditions below in particular, should be discussed with the doctor who is in charge of your treatment.
- Check the functioning of the product with your doctor or specialist retailer.
- The product is designed to treat one patient.
- Do not wear the product directly against open wounds.
- In order for the product to have a long service life and function, it may not be worn in conjunction with oily/greasy or acidic products, ointments or lotions.
- The product can come into contact with splash water.
- Please always close any Velcro fasteners present properly. Damage caused by improper closing of the Velcro fasteners does not constitute a reason for complaint.
- In the event of serious incidents when using this aid, immediately contact your doctor, your health care supply retailer or contact SPORLASTIC directly at info@sporlastic.de to report the incident.
- Any improper modifications to the product and/or improper use of the above-mentioned product exclude any product liability on the part of the manufacturer.
- Product may contain metallic objects. We recommend not putting them in contact with electromagnetic interference.
- Possible health-related mutual risks or other disadvantages with certain treatments that may arise in association with the use of the product must be discussed with the doctor in charge of treatment.

 Please return the packaging and the product to the local recycling collection system. Please comply with local regulations.

 SPORLASTIC GmbH is affiliated with a returns system for packaging and therefore fully complies with the packaging legislation in the Federal Republic of Germany.

CARE

TRIAC II can be cleaned with a damp cloth. Dry the product in the air and not in the tumble dryer, on the heater or in the microwave.

TRANSPORT AND STORAGE CONDITIONS

Please ensure that the product is stored in a dry place and is protected against moisture and sunlight. Store the product at a normal temperature and humidity.

QUALITY MANAGEMENT SYSTEM

All products from SPORLASTIC GmbH are subject to product inspection as part of our quality management system. If you still have any complaints about our product, we kindly ask that you contact your specialist retailer.

 An online version of the instructions for use can be found on our website at www.sporlastic.de.

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