

**SPORLASTIC SUPREME**  
EVEN BETTER SUPPORTS



**BETTER IN MOTION**

# SPORLASTIC HOLISTIC TREATMENT CONCEPT


## ALL-ROUND CARE

Whether it be injuries or conditions affecting the hand, elbow, spine, knee or ankle: when joints are painful, supporting and blood circulation-promoting supports in combination with accompanying training therapy measures can promote the healing process in a targeted manner.

With its SUPREME products, SPORLASTIC offers patient-centric solutions for therapy in a wide range of indications for the hand, elbow, spine, knee and ankle.


### SUPPORTS AND BRACES HELP

According to a representative **Allensbach survey**<sup>1</sup> on the effectiveness of orthopaedic aids, **76%** of wearers of supports/braces who previously had to regularly take pain-relieving medication were now able to **dispense with analgesics** either completely or partially as a result of using the aid.



Source:

1) Representative survey by the Institut für Demoskopie Allensbach on behalf of eurocom e.V. 2019 – European Manufacturers Federation for Compression Therapy and Orthopaedic Devices: “Benefits and effectiveness of medical aids”





Functional aids and appliances



Special exercise programmes



Helpful, tried and tested training equipment

## INCLUDING PHYSIOTHERAPY TO GO DIRECT BENEFIT FOR YOUR PATIENTS

SPORLASTIC SUPREME supports and braces, in combination with medical and physiotherapy measures, support a holistic treatment concept that includes the original aspects of the underlying indications. Physiotherapy exercises are part of this approach.

### Offer your patients true additional benefits:

Patient information brochures or QR-codes on selected products lead patients with smartphones or tablets to important product information as well as guided therapy exercises to promote faster recovery.

Your patients have access to this important information and exercises anytime and anywhere, and can carry out these independently, according to the motto: physiotherapy to go.



Instructed  
exercise videos

Access to important  
information via  
QR-Code anytime,  
anywhere



# SPORLASTIC SUPREME

## EVEN BETTER SUPPORTS


SPORLASTIC SUPREME products guarantee good skin compatibility thanks to carefully selected materials and tested microclimate values. The optimised, latex-free 3D flat knitted fabrics with recycled yarn adapt even better to the body's physiological movement. New comfort and stretch zones ensure even greater wearing comfort with optimum effectiveness and stabilisation.

The SPORLASTIC AIR-MATRIX silicone pads adapt perfectly to the anatomy of the bodypart. This supports the proprioceptive mode of action and an ideal fit. The perforated structure guarantees a lighter weight and better breathability.

## FOUR GOOD REASONS FOR SPORLASTIC SUPREME


- 1 Breathable 3D flat knitted fabric**  
with recycled yarn
- 2 Innovative stretch zones**  
for the perfect fit, even during movement
- 3 SPORLASTIC AIR-MATRIX silicone pads**  
breathable, lightweight, anatomically shaped
- 4 Physiotherapy to go**  
via QR CODE in every product packaging





## MALLEO-HIT® SUPREME


Ankle support with nopped  
AIR-MATRIX silicone friction pads



### INDICATIONS

- Conservative/post-operative
- Chronic capsular ligament instability
- Mild ankle sprain
- Irritation of the upper / lower ankle
- Arthrosis of the upper / lower ankle
- Rheumatological conditions
- Swelling

• **Article No.** 07174

• Colours: 

• Can be worn on the R and L

## VERTEBRADYN® SUPREME SENSO

Lumbar spine brace with  
AIR-MATRIX silicone pad for stabilisation



### INDICATIONS


- Conservative
- Lumbar pain
- Muscular imbalances
- Radicular irritation syndromes
- Spinal canal stenosis
- Preventative stabilisation

Straight cut

• **Article No.** 07545

Tailored cut


• **Article No.** 07550

• Colours: 



## GENU-HIT® SUPREME

Knee support with AIR-MATRIX silicone pad and lateral reinforcement



Soft support edge with reduced pressure

Innovative stretching zone for the perfect fit, even during movement

Lateral, stabilising spiral springs

Intarsia comfort zone for crease-free wearing


**SPORLASTIC AIR-MATRIX silicone pad:** anatomically shaped for pressure relief and guidance of the patella

Breathable 3D flat knitted fabric with recycled yarn

## INDICATIONS

- Conservative/post-operative
- Ligament injuries
- Mild ligament laxity
- Meniscal damage
- Chondropathia patellae
- Baker's cysts
- Knee arthrosis
- Post-traumatic inflammation
- Rheumatological conditions

• **Article No.** 07181

• **Colours:** 

• Can be worn on the R and L




## MANUDYN® SUPREME

Wrist support with AIR-MATRIX silicone pad for dynamic stabilisation




### INDICATIONS

- Conservative/post-operative
- Sprains
- Tendon sheath inflammation
- Wrist strain symptoms
- Following fractures
- Following removal of plaster casts
- Inflammatory conditions (post-traumatic, post-operative, rheumatism-related)

- **Article No.** 07157
- Colours: 
- Please specify R or L


## EPIDYN® SUPREME

Elbow support with nopped AIR-MATRIX silicone pads and circular strap for compression and friction



### INDICATIONS

- Conservative/post-operative
- Epicondylitis humeri radialis/ulnaris (tennis elbow, golfer's elbow)
- Arthrosis
- Ligament laxity
- Elbow sprain
- Inflammatory conditions (e.g. after fractures and operations)
- Rheumatological conditions

- **Article No.** 07337
- Colours: 
- Can be worn on the R and L

# SPORLASTIC ORTHOPAEDICS

## BETTER IN MOTION

### SIZE IN CM

#### MALLEO-HIT® SUPREME

##### Ankle circumference

1	17 - 19
2	19 - 21
3	21 - 23
4	23 - 25
5	25 - 27
6	27 - 29

#### VERTEBRADYN® SUPREME SENSO

##### Body circumference

XS	70 - 80
S	80 - 90
M	90 - 100
L	100 - 110
XL	110 - 120
XXL	120 - 130

Front height 15.5 cm  
Back height 23.0 cm

#### GENU-HIT® SUPREME

##### Measurement d f

1	25 - 28	35 - 38
2	28 - 31	38 - 41
3	31 - 34	41 - 44
4	34 - 37	44 - 47
5	37 - 40	47 - 50
6	40 - 43	50 - 53
7	43 - 46	53 - 56

#### MANUDYN® SUPREME

##### Wrist circumference

XS	14 - 15
S	15 - 16
M	16 - 17
L	17 - 18
XL	18 - 19
XXL	19 - 21

#### EPIDYN® SUPREME

##### Elbow circumference

XS	22 - 24
S	24 - 26
M	26 - 28
L	28 - 30
XL	30 - 32
XXL	32 - 34

### SPORLASTIC GERMANY

#### SPORLASTIC GmbH Medical products

Weberstrasse 1 • 72622 Nürtingen • Germany  
PO Box 14 48 • 72604 Nürtingen • Germany

Phone +49 7022 705181  
Fax +49 7022 705113  
info@sporlastic.de  
www.sporlastic.de

### SPORLASTIC QUALITY



Quality management system  
DQS-certified according to DIN  
EN ISO 13485:2016, reg. no.  
069305 MP2016



MITTELSTAND  
GLOBAL  
HEALTH MADE  
IN GERMANY

### SPORLASTIC DIGITAL



[www.facebook.com/sporlastic](http://www.facebook.com/sporlastic)



[www.youtube.com/sporlastic](http://www.youtube.com/sporlastic)



[www.instagram.com/sporlastic\\_official](http://www.instagram.com/sporlastic_official)



The product catalogue 2021 and all of our  
product brochures are also available digitally  
at [www.sporlastic.de/service/downloads](http://www.sporlastic.de/service/downloads)